

Soft, Single Patient-Use
Blood Pressure Cuffs


SunTech® *Disposable*
BP CUFFS

Providing Affordable Infection Control Safety



- + Engineered with patient comfort in mind
- + Trust in 100% leak testing

- + Assistance in standardizing your facility on one cuff connector
- + Reliability from the BP experts

For over twenty years, SunTech Medical® has focused solely on the development of blood pressure measurement technologies for the most challenging and complex clinical environments. As automated monitors have become routine for blood pressure measurement, it is important to remember that the cuff is the sensor. With this in mind, we have applied our expertise to engineer a new line of disposable blood pressure cuffs with a tapered end for enhanced patient comfort and external and internal range indicators to ensure accurate fit. This unique approach makes SunTech Medical disposable cuffs the easy choice to meet your infection control requirements at an exceptional value.



Full Range of Cuff Sizes, Including Neonatal, Infant & Child



One or Two-Tube Designs Available















+ Features & Benefits

- Patient Comfort:**
- Economical:**
- Ease of Use:**
- Patient Identification:**
- Reliable Performance:**
- Regulatory Requirements:**
- Durability:**
- Latex-Free:**
- Cuff Standardization Program (CSP):**

+ Description

- Softer materials, tapered end and rounded edges provide a gentle cuff solution.
- Priced to support single-patient use, particularly where infection control is a concern.
- Tapered end for easier cuff application. Color-coded for quick size selection. Interior & exterior range indicators ensure accurate cuff sizing and fit.
- Dedicated location for patient identification encourages single patient use for each cuff.
- 100% leak tested for trust, value, and reliability.
- Engineered to meet AHA, AAMI-SP10, and MDD guidelines for results you can trust.
- Value in a durable, single-piece design that can withstand multiple inflations.
- 100% latex-free materials promotes patient safety by minimizing risk of adverse reactions.
- Expertise in standardizing to one BP cuff connector, so your cuffs can connect to any monitor in your facility.

+ Ranges

-  Neonatal 1
3-6 cm 
-  Neonatal 2
4-8 cm 
-  Neonatal 3
6-11 cm 
-  Neonatal 4
7-13 cm 
-  Neonatal 5
8-15 cm 
-  Infant
8-13 cm 
-  Child/Child Long
12-19 cm 

- Small Adult/
Small Adult Long**
17-25 cm
- Adult/Adult Long**
23-33 cm
- Large Adult/
Large Adult Long**
31-40 cm
- Thigh**
38-50 cm

+ Symbols

-  Index Line.
-  Cuff index line must fall within range markings for an accurate measurement.
-  Artery symbol and arrow should be placed over brachial or femoral artery.
-  Product contains no latex material.
-  Product in compliance with Council Directive 93/42/EEC Medical Devices Directive.
-  Arm circumference



3046 Brecksville Road, Suite 2
Richfield, Ohio 44286

Sales & Clinical Support
800 837 5511 x200
sales@vitalcareproducts.com

www.VitalCareProducts.com



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.
Oakfield Estate
Eynsham, Oxfordshire
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235

www.SunTechMed.com

SunTech Medical, Ltd.
Level 19, Two IFC
8, Finance Street, Central
Hong Kong
Tel: +852.2251.1949
Fax: +852.2251.1950

www.CN.SunTechMed.com

SunTech Medical, Ltd.
2/F of Building A, Jinxiongda Technology Park
Guanlan, Bao'an District Shenzhen, 518110, PRC
Tel: +86-755-2958 8810
+86-755-29588986 (Sales)
+86-755-29588665 (Service)
Fax: +86-755-2958 8829

www.CN.SunTechMed.com

