


HUMIDIFIERS AND WATER TRAPS



For dependable, cost-effective performance, Salter humidifiers deliver a durable, disposable product that's ideal for use with all oxygen sources: concentrators, liquid, tanks, or wall outlets. Unique 360° rectangular diffusion ports offer quiet operation with improved hydration and less system back-pressure.

Standard Bubble Humidifiers

Available with either 3 PSI or 6 PSI pop-off alarms.

Description	Units per case	Item number
350 cc bubble humidifier 3 PSI (brown lid)	50	7100-0-50
350 cc bubble humidifier 6 PSI (black lid)	50	7600-0-50

High Flow Bubble Humidifier

Ideal for use with high-output oxygen concentrators, wall outlet source gas, or other applications where flow rates of 6 to 15 LPM are required and an audible pressure relief alert is desired. Effectively pairs with Salter 1600HF and 16SOFT-HF nasal cannulas.

Description	Units per case	Item number
High Flow 350 cc bubble humidifier (2-tone, beige and black lid)	25	7900-0-25

Humidifier Convenience Kits

Each kit includes humidifier and nasal cannula for simpler ordering and patient setup.

Description	Units per case	Item number (3 PSI)	Item number (6 to 15 LPM)
Kit with 7' soft cannula	25	—	79SOFT-7-25
Kit with 7' cannula	10	—	7907-7-10
Kit with 14' soft cannula	25	—	79SOFT-14-25
Kit with 25' cannula	20	7125-25-20	—
Kit with 25' soft cannula	10	—	79SOFT-25-10

Connector Tubing for Bubble Humidifiers

Crush resistant, 3-channel design provides confidence in uninterrupted oxygen flow. Tubing remains straight when uncoiled, allowing for greater comfort and convenience for the patient.

Description	Units per case	Item number
Concentrator humidifier adapter tubing: 1' O ₂ supply line with one molded trumpet connector end and one male thread grip end	25	9996-1-25
Concentrator humidifier adapter tubing: 15" O ₂ supply line with one molded ribbed connector end and one male thread grip end	50	SO-676
Concentrator humidifier adapter tubing: 21" O ₂ supply line with one molded ribbed connector end and one male thread grip end	50	SO-1790

Salter Inline Water Traps

A cost-effective solution for your long-term oxygen therapy patients: simple, durable construction, with design features that minimize the hazard of snagging on furniture or other obstructions.



Description	Units per case	Item number
Clear body with self-sealing end cap	25	7000-0-25
Clear body with threaded end cap	25	7001-0-25

Our most popular models are indicated by **Salter Green** throughout this catalog.