

PERCENT-O₂-LOCK® VENTURI OXYGEN SYSTEM

These masks offer an easily adjustable, variable Venturi system delivering accurate concentrations of oxygen (24%, 26%, 28%, 31%, 35%, 40%, and 50%) from a single unit. Concentration levels are easily dialed and locked in place.



Description	Units per case	Item number Elastic head strap style	Item number Elastic strap with TLC cushion	Item number Over-the-ear style
Adult air entrainment system with humidity cup, 7' O ₂ supply line, and mask	50	8150-7-50	—	8050-7-50
Adult air entrainment system with humidity cup, 7' O₂ supply line, and mask with TLC cushion	25	—	8150TLC-7-25	—
Adult trach tee air entrainment system with humidity cup and 7' O ₂ supply line	50	8054-7-50	—	—
Adult air entrainment assembly with humidity cup (without mask or O ₂ supply line tube)	50	8056-0-50	—	—
Adult air entrainment assembly without humidity cup, mask, or O ₂ supply line	50	8057-0-50	—	—
Adult air entrainment assembly with 7' O ₂ supply line and humidity cup (without mask)	50	8058-7-50	—	—
Adult air entrainment assembly with 7' O ₂ supply line (without mask or humidity cup)	50	8059-7-50	—	—
Pediatric air entrainment assembly with humidity cup, 7' O ₂ supply line, and pediatric mask	50	8121-7-50	—	—