


# ACCUTEMP®

Disposable Skin Temperature  
Probes and Probe Covers



Neonatal & Infant  
PRODUCTS

## A Cost Effective Choice

Ameritus AccuTemp probes are compatible with Air-Shields® and GE Healthcare/Ohmeda® incubators and warmers. Accurate, reliable and competitively priced; a choice you can make with confidence.

**AMERITUS**<sup>®</sup>  
Medical Products

# ACCUTEMP®

## Disposable Skin Temperature Probes and Probe Covers

### Accuracy: The Heart of Temperature Management

Clinical investigations of infants in incubators and radiant warmers clearly show the importance of accurate temperature maintenance in order to properly manage apnea, metabolic rate and insensible water loss. Ameritus AccuTemp disposable temperature probes and probe covers allow for the most accurate reading of the patient's true skin temperature.

#### Features

- **AccuTemp** disposable temperature probes are designed for single patient use during the patient's length of stay, decreasing potential risk of cross-contamination and improving hospital cost recovery
- Low profile sensor design allows probe to be held securely in place, optimizing heat conduction from the skin to the probe
- Ameritus **AccuTemp** and **AccuTemp Plus** probe covers feature highly reflective gold mylar and an added layer of foam insulation that serve as a superior barrier against both ambient and radiant heat
- **AccuTemp Plus** Insulated Hydrogel Covers provide a very gentle method of securing temperature probes to fragile skin. They may be easily and comfortably repositioned and reapplied to reduce product waste
- All temperature probes and probe covers are Latex Free and DEHP Free
- Made in the USA

#### ORDERING INFORMATION

Catalog Number	Description	Quantity
120	<b>AccuTemp</b> Acrylic Disposable Temperature Probe Cover	120/Box
120P	<b>AccuTemp</b> Acrylic Disposable Temperature Probe Cover, Preemie	
120S	<b>AccuTemp</b> Acrylic Disposable Temperature Probe Cover, Small	
175	<b>AccuTemp Plus</b> Hydrogel Disposable Temperature Probe Cover	75/Box
175-S	<b>AccuTemp Plus</b> Hydrogel Disposable Temperature Probe Cover, Small	
210	Disposable Skin Temperature Probe, Compatible with Air Shields C-100, C-200, C-450, C-550 Isolettes and PM-78 Series Warmer, IICS-90, Free Standing and Wall Mounted Birthing Room Warmers and Delivery Room Warmers	10/Box
230	Disposable Skin Temperature Probe, Compatible with GE/Healthcare/ Ohmeda® Care Plus® and Ohio® I.C. Incubators as well as Ohmeda Infant Warmer Systems 3000, 4000, and 5000 Series	10/Box
232	Disposable Skin Temperature Probe, Compatible with GE Healthcare Giraffe® Omnibed®	10/Box



**AccuTemp Disposable Temperature Probes**

Catalog Numbers: 210, 230, 232

Quantity: 10/Box



**AccuTemp Plus Insulated Hydrogel Probe Covers**

Catalog Numbers: 175-S, 175

Quantity: 75/Box



**AccuTemp Insulated Acrylic Probe Covers**

Catalog Numbers: 120, 120S, 120P

Quantity: 120/Box



Manufactured for:

Kentec Medical, Inc.

17871 Fitch, Irvine, CA 92614

TEL 949.863.0810 FAX 949.833.9730

TOLL FREE (USA only) 800.825.5996

sales@ameritusmedical.com

www.kentecmedical.com



Neonatal & Infant  
PRODUCTS

DISTRIBUTED BY:

**Vital Care Products**

**Phone:** 800-837-5511

**E-mail:** sales@vitalcareproducts.com

**Web:** vitalcareproducts.com